NTI-Kahla GmbH • Rotary Dental Instruments
 Im Camisch 3 • D-07768 Kahla/Germany
 Tel. +49(0)36 424-573-0 • Fax +49(0)36 424-573-29
 e-mail: nti@nti.de • http://www.nti.de

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GB User Information for rotating dental and surgical instruments

- Use of improper instruments can cause insufficient preparation results.
- Excessive speeds while using instruments with overly large heads endanger the stability of the instruments.
 - Excessive speeds applied near to the dental pulp can endanger its vitality. The same can be the result when using an inadequate supply of water spray.
 - Observe permissible speeds indicated on packaging.
 - The handpiece must be in proper working condition. The use of an adequate supply of water spray of a minimum of 50 ml/min on the working part is necessary in the 1,500 min⁻¹ speed range to ensure reliable cooling of the material being worked on.
 - Additional external cooling is needed for Instruments with a total length of >19 mm and a head diameter of 1.8 mm (ISO -018).
 - Moistening of the surfaces to be prepared must not be impaired by an incorrect suction extraction system. Deflection of the coolant jet can also cause heat related damage.
 - The recommended applied contact pressure range is 0.3 - 2 N (30 - 200 p) depending

- on the kind of preparation.
- Avoid tilting, levering, or blocking the instruments up by using too much contact pressure. (increased danger of fracture)
- Insert the instruments as deeply as possible without excessive use of force into the chuck and check to ensure that they are secure.
- Extremely tapered instruments will wear faster at the tip due to the smaller size.
- To maintain the longevity of these instruments, careful preparation with less contact pressure is important.
- Apply Turbo instruments over their entire grinding surface.
- Roughness obtained with a diamond coating > ISO 524 has to be finished.
- Immediately remove any instruments that are damaged, bent or no longer run concentrically.
- Depending on the application, the use of protective goggles is recommended.
- Rotating instruments have to be disinfected, cleaned, dried and sterilised before the first application on the patient and immediately after every use.

Cleaning and Disinfection				Preparation for sterilization		Sterilization	Approval for use and storage	
Instruments	Process	Substance	Exposure time					
Diamond grinding tools *validated	After mechanical cleaning, thermal disinfection and cleaning*	Water >90°C	>10 min.	Visual check and functional test; Packaging acc. to EN 868	Steam sterilization, 121°C, 15 min.* 134°C, 5 min.*	Approval or rejection in the case of an incorrect procedure to be carried out by specialist, storage in suitable sterile packaging		
TC Cutters/burs Steel burs Polishers ¹⁾ Abrasives with ceramic bond Brushes	Limitations and material compatibility must be observed After mechanical cleaning, chemical disinfection and cleaning			Disinfectant, DGHM recommended, poss. with supersonic ²⁾ concentr. acc. to manufacturer's instructions.	According to manufacturer's instructions	Remove residuals of disinfectants with water, dry. Use sterilization-resistant anticorrosives. Packaging acc. to EN 868	Steam sterilization 121°C / 15 min., 134°C / 5 min.	Approval or rejection in the case of an incorrect procedure to be carried out by specialist, storage in suitable sterile packaging

Please observe when sterilising with damp heat: DIN EN 285 and DIN EN 13060.

- Sterilisation for polishers at 121°C/ 1.1 bar for 15 min is recommended.
- The use of bur stands or trays is necessary to avoid contact damage. The use of inappropriate disinfectants and/or keeping the instruments too long in the disinfectant can result in corrosion (increased danger of breakage) and discolouration of the instruments or the colour coding can come off. No alkaline disinfectant must be used for the disinfection of silicon rubber polishers (we recommend ID 213 by Dürr Dental). Polishers must not be placed in alcohol!

Storage and Keeping Conditions: In dry conditions and protected against contaminants. Protect instruments in general against chemicals, acids, heat and extreme temperature variations. Do not use polishers after expiration date indicated.
 Disposal: In accordance with the waste directory regulation (AVV) waste directory no. 180101 or 180103 for infectious material.

Observe guidelines from Robert Koch Institut - ch. E3! (Hygiene Requirements in Dentistry)

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